



SETTING UP A PCR LABORATORY

تحضيرية تأسيس معمل بيولوجيا جزيئية

By

Prof. Dr. Asmaa Hussein


أ.د. أسماء محمد الناصر أحمد حسين

Prof. of Zoonoses & Director of the
MBRU

أستاذ الأمراض المفترضة و مدير وحدة البيولوجيا الجزيئية

Assiut University

جامعة أسيوط




The PCR laboratory should consist of **three distinct work areas**. In order to avoid the contamination problems, each area should be dedicated to a single procedure


First area  Specimen preparation occurs

Second area  Reagent preparation & PCR set-up

Third area  Amplification & detection



Created with



AREA 1: SAMPLE PREPARATION:


- ✱ **Positive-displacement pipettes or pipettors with aerosol-resistant tips**
- ✱ **Gloves & laboratory coat**
- ✱ **Refrigerator, freezer, water bath or dry -heat block laminar flow biosafety cabinet**
- ✱ **Cell lysis reagents**



Created with


AREA 2: REAGENT PREPARATION & PCR SET-UP:

- Amplification reagents & supplies
- Positive-displacement pipettes or pipettors with aerosol-resistant tips
- Laminar-flow biosafety cabinet or dead air box
- Gloves & laboratory coat
- Refrigerator & freezer
- Water bath or dry -heat block



AREA 3: AMPLIFICATION & DETECTION:

- Thermal cycler
- Pipettors with aerosol-resistant tips
- Detection equipment (electrophoresis unit, incubator, plate washer, plate reader, water bath)
- Refrigerator & freezer
- Reagents & supplies for detection




Created with

The following practices will diminish the potential for contamination:

- ★ **Each area should have dedicated supplies & reagents**
- ★ **Color coding of reagents and supplies identifies those that belong to a particular area**
- ★ **Reagents, supplies and equipment should never be taken from one area to another, three sets of pipettors are therefore essential**

★ **The workflow must be unidirectional from “clean” (pre-PCR) to “dirty” (post-PCR)**

- ★ **Dedicated labcoats and gloves should be worn at each work site; when moving to a new area, workers should put on new gloves and labcoats**




Created with

 **nitro**^{PDF} professional

download the free trial online at nitropdf.com/professional


PCR LABORATORY ORGANIZATION

Flow of samples for
PCR analysis

PRE-PCR LAB  **POST-PCR LAB**

Pre-PCR is the protocols and equipment required for the isolation of nucleic acid and the assembly of the reaction to amplify the samples

Pre-PCR Lab.




Created with

 **nitro**PDF[®] professional

download the free trial online at nitropdf.com/professional

Post-PCR Lab.



Centrifuge



Avanti Centrifuge



Optical Microscopes



Safety Cabinet

Created with



Gel tank (to place the gel mold and run electrophoresis)




Gel documentation (capture and print the image of agarose gel)




UV transilluminator (basic device for visualizing fluorescence-stained gels)




Example of an image captured by gel documentation and UV transilluminator of the gel agarose electrophoresis



Thermal Cycler



PCR microcentrifuge tubes (0.5µl)



Created with



Created with


 **nitro**PDF[®] professional

download the free trial online at nitropdf.com/professional


Strict adherence to proper laboratory technique:

- **Physically isolate PCR preparations & products**
- **Autoclave solution**
- **Aliquot reagents**
- **Use disposable gloves and change gloves often during set-up**
- **Avoid splashes**


- **Use positive-displacement pipettes or aerosol resistant tips on air-displacement pipettes**
- **“ Premix” reagents**
- **Add DNA last**
- **Choose positive & negative controls carefully**



Created with



Thank you for your attention!



Created with

 **nitro**^{PDF} professional

download the free trial online at nitropdf.com/professional